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1. Rights and responsibilities

It is the responsibility of the customer to carry out the Electrical (including the appropriate earthing), Multimedia, Neumatic, Pneumatic, Water supply and Waste installations, to the connection point of the KaVo product, using the services of a duly qualified person, according to the installation plans supplied by KaVo. The maintenance of the installation cited in the previous item is the responsibility of the owner and/or user of the product. It is not the responsibility of the authorised KaVo technician to measure or check the quality of the earthing (nor the rest of the installation). If the technician finds a fault in the earth wiring, the product will not be installed. The purchaser must pay for the travelling and overstay expenses of the technical staff involved with the installation service call outside the urban perimeter of the city where the authorised service is located. As soon as these instructions are received, the customer should contact the authorised service for clarification.

Ground characteristics
The quality of the ground structure must match the load capacity buildings according to the current standards of the site. The floor must be leveled. The chair does not require fastening to the floor. We recommend fastening the dental chair to the floor in the cases where the floor has the following features: wood, hollow blocks, depression, protuberance and unevenness.

ATTENTION!

For cases of preparation of a location where there was previously other equipment, it must be re-evaluated according to the standards requested in this document.

2. Information on the installation for the purchaser**2.1 Electrical Installation**

The distribution box has a connector where the power supply wires are connected. The purchaser must provide the electrical installation in accordance with this instructions. The electrical installation prior to the installation of the product must be done in accordance with the standards in force at the site on "Low voltage electrical installations". The wires of the power supply must have their capacity calculated taking into consideration the consumption of the product, the length of the wires and the maximum voltage drop. Wiring of the trip switch is until the product: 2.5 mm².

Maximum Consumption: 145 VA

All KaVo products must be earthed. The quality of the earthing implies the safety of the dentist, his or her auxiliary and the patient. It is recommended that the supply cables are protected by a trip switch 5A. The earthing (as well as the rest of the electrical installation which supplies the product) must be inspected, measured and corrected periodically, as it is subject to degeneration over time. The earth wire must be identified with the green yellow colouring. The neutral wire must be identified with the light blue colouring. It is recommended that the supply cables are protected by a trip switch.


Monitor Support kit

The cable signal with the video HDMI (Computer) must have female connector, 60cm out of electrical conduit. The other end (Computer), must have the output of the video signal in perfect condition. The electrical conduit transmission should not have knees 90°, but curves 90°. The KaVo assumes no liability for damage from poor quality of the installation.

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

Suction unit control wiring. Comes from where the pump is to be installed. PVC cable 2,5mm², 1m out of electrical conduit.

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


Cable signal with the video HDMI (Computer) must have female connector, 60cm out of electrical conduit.
The other end (Computer), must have female connector, 60cm out of electrical conduit.

Cable signal with the video HDMI (Computer) must have female connector, 60cm out of electrical conduit.
The other end (Computer), must have the output of the video signal.

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English

2.6 Pattern for positioning the tubes at the connection point between the equipment and the supply.
Tubes under the floor.

Referência para posicionamento do ponto de ligação entre equipamento e rede
Referencia para la posición del punto de conexiones entre el equipo y la red
Reference for positioning of connection point between equipment and the supply